


# Skin corrosion protocol

## Episkin (validated)

Model: EpiSkin™ Small Model

### DESCRIPTION



### PREDICTION MODEL

<i>In vitro</i> Result	Classification ( <i>In vivo</i> Prediction)
Mean tissue viability < 35% after 3 min exposure	<b>Corrosive:</b> Optional Sub-category 1A
Mean tissue viability ≥ 35% after 3min exposure <b>AND</b> < 35% after 60min exposure <b>OR</b> Mean tissue viability ≥ 35% after 60min exposure <b>AND</b> < 35% after 240min exposure	<b>Corrosive:</b> A combination of optional Sub-categories 1B-and-1C
Mean tissue viability ≥ 35% after 240min exposure	Non-corrosive